

geon®



10 sec

Flexible

DIGITAL

THERMOMETER

MT-F1251

Features



10 Second measurement



Memory



°C/°F Switching



Auto-off



3 in 1 Measuring ways
(Oral, Axillary, Rectal)



Fever reminder



Beep



Waterproof



Flexible tip



Low power alert



Large display



NFC wireless transmission



Specifications

Model No.	MT-F1251
Measurement range	32.0°C ~ 43.0°C (89.6°F ~ 109.4°F)
Accuracy	35.5°C ~ 42.0°C (95.9°F ~ 107.6°F) ±0.1°C (0.2°F) Other range : ±0.2°C (0.3°F)
Measuring time	10 seconds
Display unit	0.1°C or 0.1°F
Memory	Last memory
Fever reminder	Yes
Operation condition	10°C ~ 35°C , RH ≤ 85%
Storage condition	-25°C ~ 55°C , RH ≤ 85%
Battery	CR2032 battery x 1

geon® Geon Corporation

— YOU HEALTH WE CARE! —

📍 No.12, Gung Ye Road, Xihu, Chang Hua, 51448 Taiwan

✉️ Info@geon.com.tw ☎️ 886 4 8852476



Find out more

www.geonmedical.com