

**geon**<sup>®</sup>



**30 sec**

**Flexible**

**DIGITAL**

**THERMOMETER**

MT-B132FA

# Features



**30 Second measurement**



**Memory**



**°C/°F Switching**



**Auto-off**



**3 in 1 Measuring ways**  
(Oral, Axillary, Rectal)



**Fever reminder**



**Beep**



**Waterproof**



**Flexible tip**



**Low Power Alert**



## Specifications

<b>Model No.</b>	MT-B132FA
<b>Measurement range</b>	32.0°C ~ 43.0°C (89.6°F ~ 109.4°F)
<b>Accuracy</b>	35.5°C ~ 42.0°C (95.9°F ~ 107.6°F) ±0.1°C (0.2°F) Other range : ±0.2°C (0.3°F)
<b>Measuring time</b>	30 seconds
<b>Display unit</b>	0.1°C or 0.1°F
<b>Memory</b>	Last memory
<b>Fever reminder</b>	Yes
<b>Operation condition</b>	10°C ~ 35°C , RH ≤ 85%
<b>Storage condition</b>	-25°C ~ 55°C , RH ≤ 85%
<b>Battery</b>	One LR41 battery

**geon® Geon Corporation**

— YOU HEALTH WE CARE! —

 No.12, Gung Ye Road, Xihu, Chang Hua, 51448 Taiwan

 [Info@geon.com.tw](mailto:Info@geon.com.tw)  886 4 8852476



 Find out more  
[www.geonmedical.com](http://www.geonmedical.com)