

Geon®



30 sec

Rigid

DIGITAL THERMOMETER

MT-B162A

Features



30 Second measurement



Memory



°C/°F Switching



Auto-off



Beep



Fever reminder



Low Power Alert





Specifications

Model No.	MT-B162A
Measurement range	32.0°C ~ 43.0°C (89.6°F ~ 109.4°F)
Accuracy	35.5°C ~ 42.0°C (95.9°F ~ 107.6°F) ±0.1°C (0.2°F) Other range : ±0.2°C (0.3°F)
Measuring time	30 seconds
Display unit	0.1°C or 0.1°F
Memory	Last memory
Fever reminder	Yes
Operation condition	10°C ~ 35°C , RH ≤ 85%
Storage condition	-25°C ~ 55°C , RH ≤ 85%
Battery	One LR41 battery

geon® Geon Corporation

— YOU HEALTH WE CARE! —

📍 No.12, Gung Ye Road, Xihu, Chang Hua, 51448 Taiwan

✉️ Info@geon.com.tw ☎️ 886 4 8852476



Find out more

www.geonmedical.com