

**Geon**<sup>®</sup>



**30 sec**

**Digital**

**BASAL**

**THERMOMETER**

MT-B231

## Features



30 Second measurement



Memory



°C/°F Switching



Auto-off



Beep



Fever reminder



Low Power Alert



Waterproof



2 Decimal Places





## Specifications

<b>Model No.</b>	MT-B231
<b>Measurement range</b>	32.0°C ~ 43.0°C (89.6°F ~ 109.4°F)
<b>Accuracy</b>	35.5°C ~ 42.0°C (95.9°F ~ 107.6°F) ±0.1°C (0.2°F) Other range : ±0.2°C (0.3°F)
<b>Measuring time</b>	30 seconds
<b>Display unit</b>	0.1°C or 0.1°F
<b>Memory</b>	Last memory
<b>Fever reminder</b>	Yes
<b>Operation condition</b>	10°C ~ 35°C , RH ≤ 85%
<b>Storage condition</b>	-25°C ~ 55°C , RH ≤ 85%
<b>Battery</b>	One LR41 battery

**geon® Geon Corporation**

— YOU HEALTH WE CARE! —

📍 No.12, Gung Ye Road, Xihu, Chang Hua, 51448 Taiwan

✉️ [Info@geon.com.tw](mailto:Info@geon.com.tw) ☎️ 886 4 8852476



🌐 Find out more  
[www.geonmedical.com](http://www.geonmedical.com)