

Geon[®]



30 sec

Digital

BASAL

THERMOMETER

MT-B261

Features



30 Second measurement



2 Decimal Places



Memory



°C/°F Switching



Auto-off



Beep



Fever reminder



Low Power Alert



Waterproof





Specifications

Model No.	MT-B261
Measurement range	32.0°C ~ 43.0°C (89.6°F ~ 109.4°F)
Accuracy	35.5°C ~ 42.0°C (95.9°F ~ 107.6°F) ±0.1°C (0.2°F) Other range : ±0.2°C (0.3°F)
Measuring time	30 seconds
Display unit	0.1°C or 0.1°F
Memory	Last memory
Fever reminder	Yes
Operation condition	10°C ~ 35°C , RH ≤ 85%
Storage condition	-25°C ~ 55°C , RH ≤ 85%
Battery	One LR41 battery

geon® Geon Corporation

— YOU HEALTH WE CARE! —

📍 No.12, Gung Ye Road, Xihu, Chang Hua, 51448 Taiwan

✉️ Info@geon.com.tw ☎️ 886 4 8852476



🌐 Find out more
www.geonmedical.com