

Geon®



NASAL ASPIRATOR

AH2



Features



Nasal aspirator



Silicon suction head



Low noise



Low vibration



One-touch operation



Three-Stage intensity control



Easy to clean



Type C charger



Specifications

Model No.	AH2
Vacuum pressure	60kpa \pm 10%
Capacity of storage container	10 ml
Battery	3.7 V Li-ion rechargeable battery
Rated Voltage	DC 5V, 1.2A
Operation conditions	5°C ~ 40°C (41°F~104°F) / RH15% ~ 90% / 700-1060npa
Storage conditions	-25°C ~ 55°C(-13°F~158°F) / RH up to 90%
Dimension	197.2mm x 43.1mm

geon® Geon Corporation

— YOU HEALTH WE CARE! —

 No.12, Gung Ye Road, Xihu, Chang Hua, 51448 Taiwan

 Info@geon.com.tw  886 4 8852476



 [Find out more
www.geonmedical.com](http://www.geonmedical.com)