



NASAL CLEANER

H2



Features



Nasal Irrigator



Nasal aspirator



Three-Stage intensity control



Silicon suction head



Low noise



Low vibration



Easy to clean



Type C charger



Specifications

Model No.	H2
Vacuum pressure	-60kpa ± 10%
Capacity of storage container	20 ml
Spray capacity	5 ml/min
Battery	3.7 V Li-ion rechargeable battery
Operation conditions	5°C ~ 40°C (41°F~104°F) / RH15% ~ 90% / 700-1060npa
Storage conditions	-25°C ~ 70°C(-13°F~158°F) / RH up to 90%
Dimension	197.2mm x 41.3mm

geon® Geon Corporation

— YOU HEALTH WE CARE! —

 No.12, Gung Ye Road, Xihu, Chang Hua, 51448 Taiwan

 Info@geon.com.tw  886 4 8852476



Find out more

www.geonmedical.com