



NASAL ASPIRATOR



Features



Nasal aspirator



Silicon suction head



One-touch operation



Reflux prevention



Low vibration



AA battery x 2



Easy to clean



Specifications

Model No.	S2
Vacuum pressure	-60kpa ± 10%
Capacity of storage container	10 ml
Battery	2 x AA alkaline battery
Operation conditions	5°C ~ 40°C (41°F~104°F) / RH15% ~ 90% / 700-1060npa
Storage conditions	-25°C ~ 70°C(-13°F~158°F) / RH up to 90%
Dimension	36.6mm x 68.6mm x 172mm

geon® Geon Corporation

— YOU HEALTH WE CARE! —

 No.12, Gung Ye Road, Xihu, Chang Hua, 51448 Taiwan

 Info@geon.com.tw  886 4 8852476



 [Find out more
www.geonmedical.com](http://www.geonmedical.com)