

geon[®]



1 sec

Non-contact

FOREHEAD

THERMOMETER

GE-TF11



Features



1 Second measurement



°C/°F Switching



Beep



Fever reminder



Two modes switching



30 Memory Sets



Auto-off



Low Power Alert



CR2032 Battery x 1



Specifications

Model No.	GE-TF11
Measurement range	32.0°C ~ 43.0°C (89.6°F ~ 109.4°F)
Accuracy	35.0°C ~ 42.0°C ±0.2°C (95.9°F ~ 107.6°F ±0.4°F) Other range : ±0.3°C (±0.5°F)
Measuring time	1 seconds
Display unit	0.1°C or 0.1°F
Memory	30 sets
Fever reminder	Yes
Operation condition	10°C ~ 40°C , RH ≤ 85%
Storage condition	-20°C ~ 60°C , RH ≤ 85%
Battery	CR2032 battery x 1

geon® Geon Corporation

— YOU HEALTH WE CARE! —

 No.12, Gung Ye Road, Xihu, Chang Hua, 51448 Taiwan

 Info@geon.com.tw  886 4 8852476



 [Find out more
www.geonmedical.com](http://www.geonmedical.com)